




HALO[®]

OLIGO C18

Bringing the speed and efficiency of Fused-Core[®] performance to oligonucleotide separations.

A BETTER PATH TO SEPARATIONS



HALO[®] OLIGO C18

TAKING SEPARATIONS TO A HIGHER LEVEL

Built upon proven Fused-Core[®] particle technology for speed and efficiency, the HALO[®] OLIGO C18 incorporates surface modified organo-silane technology for alkaline resistance resulting in excellent stability under elevated pH operating conditions common in oligonucleotide separation methods. .


Loaded into surface passivated column hardware to address adsorption concerns, the HALO[®] OLIGO columns are ready for use with standard or bio inert instrumentation.

FEATURES OF HALO[®] OLIGO C18

- 120 Å pore size, enables separations of oligomers up to 60 bases in length
- High pH and temperature stability, designed for conditions suited for oligonucleotide separations
- UHPLC and mass spectrometry compatible stationary phase
- Surface passivated column hardware to reduce potential of stainless steel sample adsorption.

EXCELLENT LOT-TO-LOT REPRODUCIBILITY

Four different lots of HALO[®] OLIGO C18 were tested using a ladder of single stranded DNA ranging in base length from 10 mer to 60 mer.



PEAK IDENTITIES:

1. 10 mer
2. 15 mer
3. 20 mer
4. 25 mer
5. 30 mer
6. 40 mer
7. 50 mer
8. 60 mer

TEST CONDITIONS:

Columns: HALO 120 Å OLIGO C18, 2.7 µm, 2.1 x 50 mm

Part Number: P2A62-402

Mobile Phase A: 10mM TEAA, pH 8.5

Mobile Phase B: Acetonitrile

Gradient: Time	%B
0.0	5
10.0	11
11.0	11
11.5	5

Flow Rate: 0.5 mL/min

Back Pressure: 125 bar

Temperature: 60 °C

Injection Volume: 1.0 µL

Sample Solvent: 10mM Tris HCl/1mM EDTA

Detection: UV/PDA, 254 nm


Flow Cell: 1 µL

Data Rate: 100 Hz

Response Time: 0.025 sec.

LC System: Shimadzu Nexera X2


APPLICATIONS



COMPETITIVE ADVANTAGE OF HALO® OLIGO


An oligonucleotide ladder of mixed sequence and length is separated on the HALO® OLIGO C18 and a competitor oligonucleotide column under high pH conditions. The oligomers of 20 base length AND higher begin to tail significantly on the competitor column. The same oligomers show no tailing on the HALO® column demonstrating the chromatographic efficiency and speed of Fused-Core®. Note: As tailing of the competitor column could represent poor column loading, both the HALO® and competitor columns were QC tested to ensure a fair comparison.


Both columns passed QC which verified each column.



OLIGONUCLEOTIDE SEPARATION WITH LC/MS DETECTION CONDITIONS


Resolution of intermediate length synthetic oligonucleotides, and impurities, can be conducted using shallow gradients of methanol, in MS compatible IP-RP conditions. This example shows synthetic oligonucleotides of 15-21 bases resolved using 5 mM TEA as the IP reagent, with 50 mM HFIP buffer additive (left side, Abs 260 nm). In series, online ESI-MS analysis was obtained (right side), showing the charge states for the 6122 amu MW T7 promoter synthetic oligonucleotide primer. Retention of this ssDNA is determined by length and sequence (composition), which permits closely related impurities and failure sequences (n-1) to be well resolved using the HALO OLIGO C18 column.





OLIGONUCLEOTIDE PERFORMANCE MIX

By using the OLIGO C18 column under high pH conditions a sample of 6 different oligonucleotides can be separated in under 5 minutes. Using the SigmaAldrich Oligonucleotide Performance Standard Mix, the HALO® OLIGO C18 demonstrates utility as part of system suitability testing. The standard, with a range of 12 to 33 oligomers in base length, and having two 12 base length oligos, serves as an ideal performance mix.



PEAK IDENTITIES:

1. 20 mer
2. 15 mer
3. 12 mer
4. 25 mer
5. 33 mer
6. 12 mer

TEST CONDITIONS:

Column: HALO 120 Å OLIGO C18, 2.7µm, 2.1 x 50 mm

Part Number: P2A62-402

Mobile Phase A: 100mM TEAA @ pH 8.5

Mobile Phase B: Acetonitrile

Gradient:	Time	%B
	0.0	7.5
	5.0	15
	5.3	60
	5.6	60
	8.0	7.5

Flow Rate: 0.4 mL/min

Back Pressure: 142 bar

Temperature: 50 °C

Injection: 1 µL of Oligonucleotide Performance Standard Mix, 12-33 NT,

P/N: PHR8667-1EA

Sample Solvent: 10mM Tris HCl/ 1mM EDTA

Wavelength: PDA, 254 nm

Flow Cell: 1 µL


Data Rate: 100 Hz

Response Time: 0.05 sec.

LC System: Shimadzu Nexera

OLIGONUCLEOTIDE SIZE STANDARD SEPARATION

This example demonstrates the resolution of an Oligonucleotide Standard mixture (10- to 60-mer ssDNA) using triethylammonium acetate (TEAA) with absorbance detection and gradient elution with acetonitrile. This fast separation (less than 3.5 minutes) illustrates the excellent peak shape, and high resolution of oligonucleotides of up to 60 oligonucleotides in length using Advanced Material Technology's Fused-Core® technology in the 2.1 x 50 mm HALO® OLIGO C18 column.



PEAK IDENTITIES:

1. 10 mer
2. 15 mer
3. 20 mer
4. 25 mer
5. 30 mer
6. 40 mer
7. 50 mer
8. 60 mer

* Tris HCl/EDTA

TEST CONDITIONS:

Column: HALO 120 Å OLIGO C18, 2.7 µm, 2.1 x 50 mm

Part Number: P2A62-402

Mobile Phase A: 100mM TEAA, pH 8.5

Mobile Phase B: Acetonitrile

Gradient:	Time	%B
	0.0	5
	0.5	7.4
	3.5	10.7
	3.6	20
	4.1	20
	4.2	5
	9.0	5

Flow Rate: 0.5 mL/min

Back Pressure: 137 bar

Temperature: 60 °C

Injection: 2.0 µL, (10µg)

Sample Solvent: 10mM Tris HCl/1mM EDTA pH=8.0

Wavelength: PDA, 254 nm

Flow Cell: 1 µL

Data Rate: 100 Hz

Response Time: 0.025 sec.

LC System: Shimadzu Nexera X2

PRODUCT CHARACTERISTICS

Ligand: dimethyloctadecylsilane,
surface modified
Particle Size: 2.7 μm
Pore Size: 120 \AA

USP Designation: L1
Carbon Load: 5.6%
Surface Area: 75 m^2/g
Endcapped: YES

Low pH Limit: 2
High pH limit: 9
Temp limit @ low pH: 90 $^{\circ}\text{C}$
Temp limit @ high pH: 85 $^{\circ}\text{C}$

Column lifetime will vary depending on the operating temperature and the type and concentration of buffers used. Operation at extreme specifications of temperature and pH may reduce column lifetime. Consult the column Care and Use for more information.

PART NUMBERS

2.7 μm ANALYTICAL COLUMNS	
Dimensions: ID x Length (in mm)	Part Number
2.1 x 50	P2A62-402
2.1 x 100	P2A62-602
2.1 x 150	P2A62-702
4.6 x 50	P2A64-402
4.6 x 100	P2A64-602
4.6 x 150	P2A64-702

SURFACE PASSIVATED HARDWARE

INNOVATION YOU CAN TRUST – PERFORMANCE YOU CAN RELY ON

HALO®

Manufactured by:



advancedmaterialstechnology



Scantec Nordic

Analys & Mätteknik

031 336 90 00 • www.scantecnordic.se

halocolumns.com

HALO and Fused-Core are registered trademarks of Advanced Materials Technology, Inc.