


Cyclosporin A (Cyclosporin A) is a type of cyclic polypeptide antibiotic and is used as a medical drug. It is used to suppress the immune system of organ transplant patients, thereby reducing the risk of organ transplant rejection. However, due to the narrow effective therapeutic window and the potential for toxicity, it is necessary to monitor blood drug concentrations. In this study, cyclosporin A was added to human serum, extracted and purified using spin columns MonoSpin C18 and Ph, and analyzed by LC-MS/MS. (T.Kunieda)

1. Pretreatment with MonoSpin

◆ What is MonoSpin?


The MonoSpin series is a spin column using a silica monolith with uniform continuous pores. Silica monolith with high porosity is used as a carrier, so liquid can be passed through only by centrifugation. Therefore, it is possible to purify and concentrate the sample with a short and simple operation. In addition, the bed volume is small and the liquid is easily drained, so it is also suitable for small sample volumes.



how to use → <https://youtu.be/YDhZ1Dtwq3k>

application → <https://www.gls.co.jp/technique/lifescience/monospin/application/index.html>

Operating procedure



sample solution

uman serum: 0.1% formic acid: cyclosporin (50% acetonitrile containing 0.1% formic acid) = 200:190:10 μ L, centrifuged at 10,000 \times g for 1 min, and the supernatant was used as the sample solution.

* By removing fine particles by centrifugation in advance, the flow of the spin column will be improved.


Spike recovery test results

		Rec (%)	RSD(%)
0.5 ng/mL	C18	94.8	2.4
	Ph	103.8	10.0
5 ng/mL	C18	98.6	12.1
	Ph	115.6	14.8

(0.5, 5 ng/mL in sample, n=3 \times 1 days)

2.Measurement condition

Chromatogram example of standard solution




Column : InertSustain Phenyl (GL Sciences Inc.)
 (3 μm, 50 x 2.1 mm I.D.)
 Column Cat.No. : 5020-16433
 Eluent : A) CH₃CN
 B) 0.1% HCOOH in H₂O

Time (min)	A (vol%)	B (vol%)
0.0	10	90
7.0	90	10
10.0	90	10
10.1	10	90
13.0	10	90

Flow Rate : 0.4 mL/min
 Col. Temp. : 50 °C
 Detection : MS/MS(SCIEX)
 (QTRAP 6500+: ESI, Positive, SRM)
 Injection Vol. : 5 μL
 Sample : Cyclosporin A (20 ng/mL)
 Q1/Q3 = 1220/1203
 Q1/Q3 = 1220/1185

Examination of LC column and analysis conditions


■ InertSustain AQ-C18 (1.9 μm 100 x 2.1 mm I.D.)



Since the retention was strong with the ODS column and the peak was broadened, we investigated with the C8 column.




■ InertSustainSwift C8 (3 μm 100 x 3.0 mm I.D.)





Since the peak width was wide even with the C8 column, we next investigated with the phenyl column.



■ InertSustain Phenyl (3 μm 50 x 2.1 mm I.D.)



The phenyl column gave narrower peak widths, and changing the eluent gradient resulted in sharper peaks. The peak shape was further improved by increasing the column oven temperature to 50°C.



3.Related product

InertSustain Phenyl

Phenyl-bonded columns give a different separation pattern than C18 columns due to π -electron interactions. It is effective as a column of choice when separation is insufficient with a C18 column. In particular, InertSustain Phenyl has a phenyl group directly bonded to silica gel, so it has the characteristic of recognizing the difference in the electronic state of aromatic compounds to a greater extent than a general phenyl column (alkylphenyl group-bonded columns). In addition, since the phenyl group is chemically modified at high density, it has excellent stereoselectivity and exhibits unprecedented separation performance.




Column	Size	Cat.No.
InertSustain Phenyl	3 μ m, 50 x 2.1 mm I.D.	5020-16433

MonoSpin(Solid-phase extraction column for minute samples)

MonoSpin is a pretreatment spin column that uses a silica monolith with uniform continuous pores. All operations can be performed with a short centrifugation process.

The S type is ideal for samples up to 800 μ L, the L type is ideal for larger 1-8 mL samples, and the 96-well plate type is available for multi-sample processing.



MonoSpin C18

A column with octadecyl groups and reversed-phase partitioning interaction.

Ideal for extracting drugs from biological samples and desalting/concentrating peptide samples.

Product name	Quantity	Cat.No.
MonoSpin C18	50	5010-21700



MonoSpin S type

MonoSpin Ph

Reversed-phase mode type with chemically modified phenyl group. It is possible to perform pretreatment with higher selectivity by utilizing the weaker hydrophobic interaction than C18. Suitable for recovering hydrophobic drugs from biological samples.

Product name	Quantity	Cat.No.
MonoSpin Ph	50	5010-21733

GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences Inc. Japan

22-1 Nishishinjuku 6-chome
Shinjuku-ku, Tokyo
163-1130, Japan

Phone: +81-3-5323-6620
Fax: +81-3-5323-6621
Email: world@glsc.co.jp
Web: www.glsciences.com

GL Sciences Inc. USA

4733 Torrance Blvd. Suite 255
Torrance, CA 90503
USA

Phone: +1-310-265-4424
Fax: +1-310-265-4425
Email: info@glsciencesinc.com
Web: www.glsciencesinc.com

GL Sciences B.V.

Dillenburgstraat 7C
5652AM, Eindhoven
The Netherlands

Phone: +31-40-254-9531
Email: info@glsciences.eu
Web: www.glsciences.eu

GL Sciences (Shanghai) Limited

Tower B, Room 2003
Far East International Plaza
No.317 Xianxia Road, Changning District
Shanghai, China 200051

Phone: +86-21-62782272
Email: contact@glsciences.com.cn
Web: www.glsciences.com.cn



International Distributors

Visit our Website at www.glsciences.com/distributors



Analys & Mätteknik

031 336 90 00 • www.scantecnordic.se